
Simply the best fixation.



Quick release system for various surgical applications



Features:

- Interchangeable quick release head components allow various surgical applications (MIS, open surgery etc.)
- Maximum Flexibility: Just change the quick release head component according to your application
- Easy, fast and secure exchange of quick release head component with pushbutton
- Functional design of the quick release provides worry free mounting and ergonomic handling
- Quick release holder interface guarantees full compatibility for present and future quick release head component applications
- Steam sterilizable (134° C)

Simply the best fixation.



Quick release holder

without quick release adapter and rails clamp



Order no.	Description	Head component Order no.	Articulated arm Order no. Length	Column Order no. Length	Weight
3840.95	Quick release holder (QR), straight column	5.700	2.760 400 mm	5.468-05 400 mm	1.75 kg
3840.96	Quick release holder (QR), L-column	5.700	2.760 400 mm	5.495-05 450×100 mm	2.05 kg

Subject to modifications

Quick release head component

Designed for various surgical applications





5.7

Order no.	Description	Application	Connection	Range	Weight kg
5.716	Varioholder (QR)	Round instruments	Quick release adapter	Ø 4.5–12.5 mm	0.18
5.731	Combined instrument holder (QR)	Instruments	Quick release adapter	Ø 4–18 mm Opening 0–18 mm	0.19

QR = Quick Release Subject to modifications

Rails clamp





!

Order no.	Description	Material	Rail dimension	Column	Weight
5.752	Radial rails clamp	Stainless steel	Height: 25-32 mm Width: 9-10 mm	For revolving and height-adjustable column Ø 16 mm	0.90 kg
5.780	Rails clamp	Stainless steel	Height: 25-32 mm Width: 9-10 mm	For height-adjustable column Ø 16 mm	0.90 kg

Subject to modifications

Baitella AG

Thurgauerstrasse 70, CH-8050 Zurich, Switzerland Phone +41 44 305 80 00, Fax +41 44 305 80 05 info@baitella.com, www.fisso.com



